



Recombinant CA125 (MUC16) (DAG-WT1833)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90 % as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, LFIA,
Format	Liquid or Lyophilized powder
Size	100 µg, 500 µg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C. Avoid repeated freeze/thaw cycles

BACKGROUND

Introduction	CA-125 has found application as a tumor marker or biomarker that may be elevated in the blood of some patients with specific types of cancers, or other benign conditions. MUC16 (CA125) is a serum marker that is used routinely in gynecologic practice to monitor patients with ovarian cancer.
Keywords	CA125; MUC16; Tumor marker