



Recombinan Epididymis Protein 4 (HE4) (DAG-WT184)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 85% (SDS-PAGE)
Conjugate	His
Applications	Immunoassays
Molecular Weight	10,036
Format	Liquid
Concentration	> 0.5 mg/mL
Size	100 µg, 1 mg, 10 mg
Buffer	Solution in 20 mM Tris, 0.5 M sodium chloride, pH 8.0, with 50% glycerol (v/v)
Preservative	None
Storage	Store at -20°C for long term. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Human epididymis protein 4 (HE4) belongs to the family of whey acidic four-disulfide core (WFDC) proteins. Currently, the biologic function of HE4 is unknown.
Keywords	EDDM4 protein;HE-4 protein;HE 4 protein;WAP5 protein;WFDC2 protein