



Recombinant Human LH alpha/beta heterodimer protein (DAG-WT1844)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90 % as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, LFIA,
Format	Liquid or Lyophilized powder
Size	100 µg, 500 µg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C. Avoid repeated freeze/thaw cycles

BACKGROUND

Introduction	Luteinizing hormone is a hormone produced by the anterior pituitary gland. In females, an acute rise of LH called the LH surge triggers ovulation and development of the corpus luteum. In males, where LH had also been called interstitial cell-stimulating hormone (ICSH), it stimulates Leydig cell production of testosterone.
Keywords	Luteinizing hormone; LH