



# Native Hemopexin (HPX) (DAG-WT185)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 95% (SDS-PAGE)
Conjugate	Unconjugated
Cas.No	9013-71-2
Applications	Immunoassays
Molecular Weight	57,000
Format	Lyophilized
Size	1 mg, 10 mg, 100 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C for long term. Avoid freeze / thaw cycles.

## BACKGROUND

Introduction	Hemopexin (or haemopexin; Hpx; Hx), also known as beta-1B-glycoprotein, is a glycoprotein that in humans is encoded by the HPX gene and belongs to hemopexin family of proteins.
Keywords	Haemopexin; HPX; beta-1B-glycoprotein

## GENE INFORMATION

Gene Name	HPX
-----------	-----