



# Recombinant Human NGAL protein (DAG-WT1850)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Purity</b>	> 90 % as determined by SDS-PAGE
<b>Conjugate</b>	His
<b>Applications</b>	ELISA, CLIA, LFIA,
<b>Molecular Weight</b>	24 kDa
<b>Format</b>	Liquid or Lyophilized powder
<b>Size</b>	100 µg, 500 µg
<b>Buffer</b>	PBS, pH 7.4
<b>Preservative</b>	None
<b>Storage</b>	Store at -20°C to -80°C. Avoid repeated freeze/thaw cycles

## BACKGROUND

<b>Introduction</b>	LCN2 (lipocalin 2) is a protein-coding gene. Diseases associated with LCN2 include appendicitis, and proctitis. GO annotations related to this gene include iron ion binding and protein homodimerization activity. An important paralog of this gene is LCN12.
<b>Keywords</b>	LCN2; Lipocalin 2; NGAL