



Recombinant EBV Rta [His] (DAG-WT1853)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90 % as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, LFIA,
Format	Liquid
Size	1 mg
Buffer	25mM Tris-HCl, 0.02%SDS, pH 8.5
Preservative	None
Storage	Store at -20°C to -80°C. Avoid repeated freeze/thaw cycles

BACKGROUND

Introduction	The Epstein-Barr virus (EBV), also called Human herpes virus 4 (HHV-4), is a virus of the herpes family (which includes Herpes simplex virus and Cytomegalovirus). On infecting the B-lymphocyte, the linear virus genome circularizes and the virus subsequently persists within the cell as an episome. The virus can execute several distinct programs of gene expression which can be broadly categorized as being lytic cycle or latent cycle.
Keywords	Epstein-Barr virus; EBV; Rta