



Tropicamide [BSA] (DAG-WT1857)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA,
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	Tropicamide is used to dilate (enlarge) the pupil so that the doctor can see into the back of your eye. It is used before eye examinations, such as cycloplegic refraction and examination of the fundus of the eye. Tropicamide may also be used before and after eye surgery.
Keywords	Tropicamide; Mydriacyl