



Prometryn [BSA] (DAG-WT1865B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA,
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	Prometryn is a diamino-1,3,5-triazine that is N,N'-di(propan-2-yl)-1,3,5-triazine-2,4-diamine substituted by a methylsulfaneyl group at position 6. It has a role as a herbicide, a xenobiotic and an environmental contaminant. It is a diamino-1,3,5-triazine and a methylthio-1,3,5-triazine.
Keywords	Prometryn; 7287-19-6; Prometryne; Gesagard; Prometrex