



Isoprothiolane [OVA] (DAG-WT1874O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA,
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	Isoprothiolane is a highly systemic fungicide and IGR that inhibits the phospholipid biosynthesis (methyltransferase/group 2).
Keywords	Isoprothiolane; Fungicide