



Profenofos [BSA] (DAG-WT1886B)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA, LFIA,
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	An insecticide/miticide, profenofos is a restricted use pesticide used solely on cotton, to control the tobacco budworm, cotton bollworm, armyworm, cotton aphid, whiteflies, spider mites, plant bugs, and fleahoppers.
Keywords	Profenofos; Insecticide