



## Cefoperazone [BSA] (DAG-WT1891B)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                  |
|--------------|------------------|
| Conjugate    | BSA              |
| Applications | ELISA, LFIA,     |
| Format       | Liquid           |
| Size         | 1 mg             |
| Buffer       | 10mM PBS, pH 7.4 |
| Preservative | None             |
| Storage      | Store at -20°C   |
| Ship         | Wet ice          |

### BACKGROUND

|              |                                                                                                                                                                                                                                                        |
|--------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | Cefoperazone injection is used to treat bacterial infections in the different parts of the body. Cefoperazone injection belongs to the class of medicines known as cephalosporin antibiotics. It works by killing bacteria or preventing their growth. |
| Keywords     | Cefoperazone; Cephalosporin                                                                                                                                                                                                                            |