



Fenfluramine [OVA] (DAG-WT1902O)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	OVA
Applications	ELISA, LFIA,
Format	Liquid
Size	1 mg
Buffer	10mM PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	Fenfluramine is fenfluramine is a phenethylamine that is structurally similar to serotonin. Due to its ability to increase extracellular serotonin levels, modulate serotonergic and other neurologic receptors, and control neurotransmission, it is effective in treating pharmacoresistant seizures.
Keywords	Fenfluramine; Fintepla