



## Methylone [BSA] (DAG-WT192)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Product Overview	Methylone conjugated with BSA
Conjugate	BSA
Applications	LFIA/ELISA
Format	Liquid
Size	1 mg
Buffer	0.01M pH7.2 PBS
Preservative	0.01% Sodium Azide
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.

### BACKGROUND

Introduction	Methylone is an empathogen and stimulant psychoactive drug. It is a member of the substituted amphetamine, substituted cathinone and substituted methylenedioxyphenethylamine classes.
Keywords	Methylone; MDMA; Psychoactive drug; Drugs of Abuse