



Recombinant Human Osteopontin [His] (aa 1-314) (DAG-WT1922)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	NP_001035147.1
Purity	> 80% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Molecular Weight	35 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	Osteopontin (OPN) is one of the noncollagenous proteins present in bone matrix. The presence of OPN in various organs, including those with and without matrix, and also in plasma, suggests that this molecule could act both as a structural molecule and as a humoral factor, or cytokine.
Keywords	Osteopontin; OPN