



Recombinant Human β -catenin [His, GST] (DAG-WT1925)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	P35222
Purity	> 85% as determined by SDS-PAGE
Conjugate	GST
Applications	WB, ELISA
Molecular Weight	116 kDa
Format	Liquid
Size	100 μ g
Buffer	50mM Tris, 150mM NaCl, 25% glycerol, pH 8.0, 0.1mM EDTA, 1mM TCEP, 0.4mM PMSF, 0.5mM GSH
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Beta-catenin (β -catenin) is a multifunctional, 90 kD protein that contributes to cell development under normal physiological conditions [55]. β -Catenin is a crucial transcriptional factor in Wingless–Int (Wnt) signaling, and plays important role in stem cell renewal and organ regeneration
Keywords	Beta-catenin; β -catenin; Catenin beta-1