



Recombinant Human Sclerostin Protein [His] (DAG-WT1932)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	NP_079513.1
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Molecular Weight	30 kDa
Format	Lyophilized
Size	50 µg, 500 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C
Ship	Wet ice

BACKGROUND

Introduction	Sclerostin is a protein containing 188 amino acids with two N-glycosylation sites. It is secreted by osteocytes and inhibits the genesis of osteoblasts, thereby resulting in a decrease in bone formation. Circulating sclerostin levels are higher in postmenopausal than in premenopausal women.
Keywords	Sclerostin; SOST