



Risperidone [BSA] (DAG-WT194)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	Risperidone conjugated with BSA
Conjugate	BSA
Applications	LFIA/ELISA
Format	Liquid
Size	1 mg
Buffer	0.01M pH7.2 PBS
Preservative	0.01% Sodium Azide
Storage	Store at -20°C for long term. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Risperidone is a second-generation antipsychotic (SGA) medication used in the treatment of a number of mood and mental health conditions including schizophrenia and bipolar disorder. It is one of the most widely used SGAs. Paliperidone, another commonly used SGA, is the primary active metabolite of risperidone (i.e. 9-hydroxyrisperidone).
Keywords	Risperidone; Perseris; Risperdal; Antipsychotic