



Native Vitamin D binding protein (DAG-WT1944)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95%
Conjugate	N/A
Applications	WB, ELISA
Molecular Weight	56 kDa
Format	Lyophilized
Size	1 mg
Buffer	Lyophilized from 20 mM Na phosphate, pH 7.4, with 150 mM NaCl
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Vitamin D-binding protein (DBP), also/originally known as gc-globulin (group-specific component), is a protein that in humans is encoded by the GC gene. DBP is genetically the oldest member of the albuminoid family and appeared early in the evolution of vertebrates.
Keywords	Vitamin D-binding protein; DBP