



Recombinant Human CALM2 Protein [His] (DAG-WT1945)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	P0DP24
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	WB, ELISA
Molecular Weight	22 kDa
Format	Lyophilized
Size	100 µg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	This gene is a member of the calmodulin gene family. There are three distinct calmodulin genes dispersed throughout the genome that encode the identical protein, but differ at the nucleotide level. Calmodulin is a calcium binding protein that plays a role in signaling pathways, cell cycle progression and proliferation.
Keywords	CALM2; Calmodulin 2