



# Recombinant Human CALCB Protein [GST] (DAG-WT1947)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	P10092
Purity	> 90% as determined by SDS-PAGE
Conjugate	GST
Applications	WB, ELISA
Molecular Weight	40.7 kDa
Format	Liquid or Lyophilized powder
Size	20 µg, 100 µg
Buffer	Tris/PBS-based buffer
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Introduction	Calcitonin gene-related peptide (CGRP) is a member of the calcitonin family of peptides consisting of calcitonin, amylin, adrenomedullin, adrenomedullin 2 (intermedin) and calcitonin -receptor-stimulating peptide. Calcitonin is mainly produced by thyroid C cells whilst CGRP is secreted and stored in the nervous system.
Keywords	Calcitonin gene-related peptide 2; CALCB