



Inactivated SARS-CoV-2 Wild-Type (Hubei) Culture Fluid (DAG-WT1975)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^5$ TCID ₅₀ /mL
Format	Liquid
Size	100 μ L
Buffer	PBS buffer
Preservative	None
Storage	Store at 2°C to 8°C

BACKGROUND

Introduction	SARS-CoV-2 is the virus responsible for COVID-19. It was previously known as 2019 Novel coronavirus (2019-nCoV), is a positive-sense single-stranded RNA virus. SARS-CoV-2 has close genetic similarity to bat coronaviruses, from which it likely originated. An intermediate animal reservoir such as a pangolin is also thought to be involved in its introduction to humans. From a taxonomic perspective, SARS-CoV-2 is classified as a strain of the species severe acute respiratory syndrome-related coronavirus (SARSr-CoV).
Keywords	SARS-CoV-2; Coronavirus; COVID-19; Wild-Type