



Inactivated *S. epidermidis* Culture Fluid (DAG-WT1979)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Control
Bio-activity	$\geq 10^6$ CFU/mL
Format	Liquid
Size	1 mL
Buffer	DMEM
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Staphylococcus epidermidis is a Gram-positive bacterium, and one of over 40 species belonging to the genus Staphylococcus. It is part of the normal human microbiota, typically the skin microbiota, and less commonly the mucosal microbiota and also found in marine sponges.
Keywords	Staphylococcus epidermidis; <i>S. epidermidis</i>