



Recombinant Human PS20 Protein [His] (DAG-WT1987)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	25 kDa
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords PS20; WFDC1; WAP four-disulfide core domain protein 1; Prostate stromal protein ps20

GENE INFORMATION

UniProt ID Q9HC57