



Recombinant Human Alkaline Phosphatase Protein (Placental type) [His] (DAG-WT1988)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	53.8 kDa
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords Alkaline Phosphatase; ALP; ALPP; Alkaline phosphatase Regan isozyme; PLAP-1

GENE INFORMATION

UniProt ID P05187