



Biotinylated Recombinant Human PLAU/uPA Protein (active form) [His, Avi] (DAG-WT2002)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His, Avi
Applications	Immunoassays
Molecular Weight	57.6 kDa
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords Urokinase-type plasminogen activator; PLAU

GENE INFORMATION

UniProt ID P00749