



# Recombinant Human Hepcidin/HAMP Protein [GST] (DAG-WT2004)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	GST
Applications	Immunoassays
Molecular Weight	29.08 kDa
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

## BACKGROUND

**Keywords** Hepcidin; HAMP; LEAP-1; PLTR

## GENE INFORMATION

**UniProt ID** P81172