



# Recombinant Human PGLYRP1 Protein [hFc] (DAG-WT2007)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                  |  |
|------------------|--|
| Species          | Human  |
| Purity           | > 95% as determined by SDS-PAGE  |
| Conjugate        | hFc  |
| Applications     | Immunoassays   |
| Molecular Weight | 46.2 kD  |
| Format           | Lyophilized  |
| Size             | 100 µg, 500 µg   |
| Buffer           | Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization |
| Preservative     | None   |
| Storage          | Store at -20°C   |

## BACKGROUND

**Keywords** Peptidoglycan recognition protein 1; PGLYRP1; TNFSF3L

## GENE INFORMATION

**UniProt ID** O75594