



# Recombinant Human NGAL/Lipocalin-2 Protein [hFc] (DAG-WT2021)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	hFc
Applications	Immunoassays
Molecular Weight	47.2 kD
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Introduction	LCN2 (lipocalin 2) is a protein-coding gene. Diseases associated with LCN2 include appendicitis, and proctitis. GO annotations related to this gene include iron ion binding and protein homodimerization activity. An important paralog of this gene is LCN12.
Keywords	LCN2; Lipocalin 2; NGAL

# GENE INFORMATION

UniProt ID	P80188-1
------------	----------

---