



Biotinylated Recombinant Human IFN gamma/IFNG Protein [His, Avi] (DAG-WT2039)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His, Avi
Applications	Immunoassays
Molecular Weight	19 kDa
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Interferon-gamma (IFN-γ) is a cytokine that plays an important role in inducing and modulating an array of immune responses.
Keywords	Interferon gamma; IFNG

GENE INFORMATION

