



Recombinant allergen Tri a 19 [His] (DAG-WT211)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Tri a 19, or omega-5 gliadin, the major allergen of wheat, has been prepared as a recombinant mature allergen fused to a his-tag.
Purity	> 95 % as determined by SDS-PAGE.
Conjugate	His
Applications	Immunoassays
Molecular Weight	66 kDa
Format	Liquid
Size	1 mg
Buffer	20 mM phosphate buffer pH 6.0 and 0.15 M NaCl, 1 M urea and 0.1% polyoxyethylene (10) tridecyl ether
Preservative	None
Storage	Store at -20°C to -80°C. Avoid multiple freeze/thaw cycles.

BACKGROUND

Keywords Phl p 5a; Phleum pratense; Grass pollen