



Biotinylated Human CA125 Protein [His, Avi] (DAG-WT2135)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His, Avi
Applications	Immunoassays
Molecular Weight	31.3 kDa
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	CA-125 has found application as a tumor marker or biomarker that may be elevated in the blood of some patients with specific types of cancers, or other benign conditions. MUC16 (CA125) is a serum marker that is used routinely in gynecologic practice to monitor patients with ovarian cancer.
Keywords	CA125; MUC16; Tumor marker

GENE INFORMATION

UniProt ID	Q8WXI7
------------	--------
