



Biotinylated Recombinant Human FGFR1 alpha Protein [His, Avi] (DAG-WT2150)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His, Avi
Applications	Immunoassays
Molecular Weight	42 kDa
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Tyrosine-protein kinase that acts as a cell-surface receptor for fibroblast growth factors and plays an essential role in the regulation of embryonic development, cell proliferation, differentiation and migration. Required for normal mesoderm patterning and correct axial organization during embryonic development, normal skeletogenesis and normal development of the gonadotropin-releasing hormone (GnRH) neuronal system.
---------------------	---

Keywords

Fibroblast growth factor receptor 1; FGFR1; FLT-2; CD331

GENE INFORMATION

UniProt ID

P11362
