



## Recombinant Asp f 5 [His] (DAG-WT2178)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                                 |
|--------------|---------------------------------|
| Species      | Aspergillus fumigatus           |
| Purity       | > 90% as determined by SDS-PAGE |
| Conjugate    | His                             |
| Applications | ELISA                           |
| Format       | Liquid                          |
| Size         | 1 mg                            |
| Buffer       | PBS buffer                      |
| Preservative | None                            |
| Storage      | Store at -20°C                  |

### BACKGROUND

|          |                                |
|----------|--------------------------------|
| Keywords | Asp f 5; Aspergillus fumigatus |
|----------|--------------------------------|