



## Recombinant Lit v 3 [His] (DAG-WT2186)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	Litopenaeus vannamei
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Store at -20°C

### BACKGROUND

Keywords	Lit v 3; Litopenaeus vannamei
----------	-------------------------------