



Recombinant Sola I 4 [His] (DAG-WT2193)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Solanum lycopersicum
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Sola I 4; Solanum lycopersicum
----------	--------------------------------