



Recombinant Human ApoA2 [His] (DAG-WT2206)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	9.74 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	May stabilize HDL (high density lipoprotein) structure by its association with lipids, and affect the HDL metabolism
Keywords	Apolipoprotein A-II; APOA2

GENE INFORMATION

