



Recombinant Human PG I [His] (DAG-WT2208)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	41.3 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent
Keywords	Pepsin A-5; PGA5

GENE INFORMATION

UniProt ID	P0DJD9
------------	--------