



Recombinant Human GOLM1 [His] (DAG-WT2217)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	42.6 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM PB, 150mM NaCl, pH7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Golgi membrane protein 1 (GOLM1) also known as Golgi phosphoprotein 2 or Golgi membrane protein GP73 is a protein that in humans is encoded by the GOLM1 gene. Two alternatively spliced transcript variants encoding the same protein have been described for this gene.
Keywords	Golgi membrane protein 1; GOLM1

GENE INFORMATION

