



Recombinant Human GDF15-HSA [His] (DAG-WT2222)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	81 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	10mM NaAc, pH 4.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Brainstem-restricted receptor for GDF15 which regulates food intake, energy expenditure and body weight in response to metabolic and toxin-induced stresses.
Keywords	Growth differentiation factor 15; GDF15; Cardiovascular

GENE INFORMATION

