



Recombinant Human sTR [His] (DAG-WT2225)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	75 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM PB, 150mM NaCl, 5% Trehalose, 5% Mannitol, 0.01% Tween80, pH7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Transferrin receptor 1 is a cellular receptor for human heme-albumin. Protein Expression of Cyclin B1, Transferrin Receptor, and Fibronectin Is Correlated with the Prognosis of Adrenal Cortical Carcinoma.
Keywords	Transferrin receptor protein 1; TR

GENE INFORMATION

