



Recombinant Human Serpin E1 [His] (DAG-WT2226)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	43.8 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM HAc-NaAc, 150mM NaCl, pH 4.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Plasminogen activator inhibitor-1 (PAI-1) also known as endothelial plasminogen activator inhibitor (serpin E1) is a protein that in humans is encoded by the SERPINE1 gene. Elevated PAI-1 is a risk factor for thrombosis and atherosclerosis.
Keywords	Plasminogen activator inhibitor 1; SERPINE1

GENE INFORMATION

