



Recombinant Human TPO [His] (DAG-WT2227)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	92.8 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM Tris-HCl, 500mM NaCl, 10% Glycerol, pH 8.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Thyroid peroxidase (TPO) is an enzyme normally found in the thyroid gland. TPO plays an important role in the production of thyroid hormones. A TPO test detects antibodies against TPO in the blood.
Keywords	Thyroid peroxidase; TPO

GENE INFORMATION

