



Recombinant Human IFN γ (DAG-WT2230)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	16.9 kDa
Format	Liquid
Size	50 μ g, 500 μ g, 1 mg
Buffer	20mM PB, 300mM NaCl, 5% Mannitol, 0.1% Tween80, 5% Trehalose, pH6.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Interferon-gamma (IFN- γ) is a cytokine that plays an important role in inducing and modulating an array of immune responses.
Keywords	Interferon gamma; IFNG

GENE INFORMATION

UniProt ID	P01579
------------	--------