



Recombinant Human CRP Protein [His] (DAG-WT2232)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	24 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM Tris-HCl, 280mM NaCl, 5mM CaCl ₂ , pH8.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	<p>C-reactive protein (CRP) is a protein found in the blood, the levels of which rise in response to inflammation (i.e. C-reactive protein is an acute-phase protein). Its physiological role is to bind to phosphocholine expressed on the surface of dead or dying cells (and some types of bacteria) in order to activate the complement system via the C1q complex. CRP is synthesized by the liver in response to factors released by fat cells (adipocytes). It is a member of the pentraxin family of proteins. It is not related to C-peptide or protein C. C-reactive protein was the first pattern recognition receptor (PRR) to be identified.</p>
---------------------	---

KeywordsCRP; C-reactive protein

GENE INFORMATION**UniProt ID**P02741
