



# Recombinant Human Haptoglobin [His] (DAG-WT2234)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|                  |                                 |
|------------------|---------------------------------|
| Species          | Human                           |
| Purity           | > 95% as determined by SDS-PAGE |
| Conjugate        | His                             |
| Applications     | Immunoassays                    |
| Molecular Weight | 15.9&28.3kDa                    |
| Format           | Liquid                          |
| Size             | 50 µg, 500 µg, 1 mg             |
| Buffer           | PBS, pH 7.4                     |
| Preservative     | None                            |
| Storage          | Store at -20°C                  |

## BACKGROUND

|              |   |
|--------------|---|
| Introduction | <p>Haptoglobin is a protein made by your liver. It attaches to a certain type of hemoglobin. Hemoglobin is a protein in your red blood cells that carries oxygen from your lungs to the rest of your body. Most hemoglobin is located inside red blood cells, but small amounts circulate in the bloodstream.</p> |
| Keywords     | Haptoglobin; HP   |

# GENE INFORMATION

UniProt ID

P00738

---