



Recombinant Human NT-3 (DAG-WT2235)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	13.62 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM PB, 250mM NaCl, pH 7.2
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Neurotrophin-3 (NT-3) is a protein that plays a vital role in the survival and growth of nerve cells in the brain and other parts of the nervous system. It belongs to the family of neurotrophins, which are molecules that support the development and maintenance of neurons.
Keywords	Neurotrophin-3; NTF3

GENE INFORMATION

