



Recombinant Human SPP1 [His] (DAG-WT2236)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	34.7 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Osteopontin (OPN) is one of the noncollagenous proteins present in bone matrix. The presence of OPN in various organs, including those with and without matrix, and also in plasma, suggests that this molecule could act both as a structural molecule and as a humoral factor, or cytokine.
Keywords	Osteopontin; OPN

GENE INFORMATION

