



## Recombinant Human Jo-1 (DAG-WT2237)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	57.5 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

### BACKGROUND

Introduction	Anti-Jo-1 autoantibodies which recognize histidyl-tRNA synthetase identify patients with the rare rheumatologic disease, anti-histidyl-tRNA synthetase syndrome (Jo-1 ARS), a phenotypically distinct subset of idiopathic inflammatory myopathies (IIM)
Keywords	Histidine--tRNA ligase; HARS1; Jo-1

### GENE INFORMATION

