



Recombinant Dermatophagoides pteronyssinus Der p 2 (DAG-WT2242)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	14.9 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	PBS, 50% Glycerol, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	Dermatophagoides pteronyssinus; Der p 2
----------	---

GENE INFORMATION

UniProt ID	P49278
------------	--------