



Recombinant Human ApoH [His] (DAG-WT2243)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	37.3 kDa
Format	Liquid
Size	50 µg, 500 µg, 1 mg
Buffer	20mM PB, 150mM NaCl, pH 7.2
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Binds to various kinds of negatively charged substances such as heparin, phospholipids, and dextran sulfate. May prevent activation of the intrinsic blood coagulation cascade by binding to phospholipids on the surface of damaged cells.
Keywords	ApoH; Beta-2-glycoprotein 1

GENE INFORMATION

